



Special Applications with **KENTAMED hARTT 400**




CE 0044


The electrocautery unit **hARTT 400** is the newest representative of the **hARTT family**, designed and produced by **KENTAMED Ltd.** Besides its raised maximal power, the new **hARTT 400** puts special accent on bipolar cutting and coagulation in wet environment (SALINE).

Bipolar Transurethral Resection (TUR)

In combination with a bipolar resectoscope, **KENTAMED hARTT 400** is a perfect solution for transurethral resection of the prostate as well as resection of bladder tumors.




Bipolar Resectoscope




Bipolar Resection of the prostate



Bipolar Coagulation of the prostate




Bipolar Resection of bladder tumors



Arthroscopy

KENTAMED hARTT 400 is the perfect choice for arthroscopic surgery. It certainly provides excellent performance when used for monopolar arthroscopy but the bipolar application is where this model really shines. With the special set of accessories provided by our partners, **hARTT 400** can execute arthroscopic textbook cutting, coagulation and vaporization in saline.

Bipolar Arthroscopy in Saline



Monopolar Arthroscopy



Manufacturer:

KENTAMED Ltd.

41, Nestor Abadjiev St.

POB 203, Plovdiv 4023, Bulgaria

tel./fax: +359 32 622 934

e-mail: sales@kentamed.com

www.kentamed.com



Distributor:

